


# Merging Cell Biology & Technology

## Mem-Lok<sup>®</sup> Resorbable Collagen Membrane



### **Easy Handling, Long Lasting:**

- Predictable resorption rate of 26-38 weeks
- Cell occlusive membrane prevents soft tissue ingrowth
- Macromolecular pore size allows vital nutrient transfer
- Favorable handling characteristics adapt to various deficiencies
- Highly purified Type I Collagen derived from bovine achilles tendon

### **Ordering Information**

- RCM-ML1520 Resorbable Collagen Membrane 15mm x 20mm
- RCM-ML2030 Resorbable Collagen Membrane 20mm x 30mm
- RCM-ML3040 Resorbable Collagen Membrane 30mm x 40mm

For more information or to place an order, call  
1-888-246-8338 or contact your local representative.

Mem-Lok is manufactured by Collagen Matrix, Inc.  
Mem-Lok is not available in all markets.  
References available upon request.

**BIOHORIZONS<sup>®</sup>**  
mem-lok